

MARK3 Rabbit mAb

Catalog No: #56483

Package Size: #56483-1 50ul #56483-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MARK3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Mouse, Human
Specificity	MARK3 Antibody detects endogenous levels of total MARK3
Immunogen Description	A synthesized peptide derived from human MARK3
Other Names	C-TAK1; cTAK1; Emk2; KP78; Mark3; PAR1A;
Accession No.	Uniprot:P27448
Uniprot	P27448
Calculated MW	86kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000ICC/IF:1:50~1:200IP:1:30FC:1:50

Images



Western blot analysis of MARK3 expression in K562 cell lysate.

Product Description

Involved in the specific phosphorylation of microtubule-associated proteins for tau, MAP2 and MAP4. Phosphorylates CDC25C on 'Ser-216'.

Note: This product is for in vitro research use only