

IRF-9 Rabbit mAb

Catalog No: #56496

Package Size: #56496-1 50ul #56496-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	IRF-9 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human
Specificity	IRF-9 Antibody detects endogenous levels of total IRF-9
Immunogen Description	A synthesized peptide derived from human IRF-9
Other Names	Irf9; ISGF3; ISGF3G; p48;
Accession No.	Uniprot:Q00978
Uniprot	Q00978
Calculated MW	48kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:1000~1:5000|ICC/IF:1:50~1:200|FC:1:50

Images



Western blot analysis of Interferon regulatory factor 9 expression in Jurkat cell lysate.

Product Description

Transcription regulatory factor that mediates signaling by type I IFNs (IFN-alpha and IFN-beta). Following type I IFN binding to cell surface receptors, Jak kinases (TYK2 and JAK1) are activated, leading to tyrosine phosphorylation of STAT1 and STAT2.

Note: This product is for in vitro research use only