

GOLPH4 Rabbit mAb

Catalog No: #56511

Package Size: #56511-1 50ul #56511-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	GOLPH4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	GOLPH4 Antibody detects endogenous levels of total GOLPH4
Immunogen Description	A synthesized peptide derived from human GOLPH4
Other Names	GIMPC; golim4; GOLPH 4; GPP130; P 138; P138;
Accession No.	Uniprot:O00461
Uniprot	O00461
Calculated MW	130kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000


IHC: 1:100~1:500

Images



Western blot analysis of GOLPH4 expression in HepG2 cell lysate.

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Product Description

Plays a role in endosome to Golgi protein trafficking; mediates protein transport along the late endosome-bypass pathway from the early endosome to the Golgi.

Note: This product is for in vitro research use only