

EBF1 Rabbit mAb

Catalog No: #56515

Package Size: #56515-1 50ul #56515-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	EBF1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse
Specificity	EBF1 Antibody detects endogenous levels of total EBF1
Immunogen Description	A synthesized peptide derived from human EBF1
Other Names	COE1; EBF; Ebf1; O/E 1; OE1; OLF1;
Accession No.	Uniprot:Q9UH73
Uniprot	Q9UH73
Calculated MW	64kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:1000

Images



Western blot analysis of EBF1 expression in Raji cell lysate.

Product Description

Transcriptional activator which recognizes variations of the palindromic sequence 5'-ATTCCNNGGAATT-3'.

Note: This product is for in vitro research use only