

TNNI2 Rabbit mAb

Catalog No: #56516

Package Size: #56516-1 50ul #56516-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TNNI2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human Mouse Rat
Specificity	TNNI2 Antibody detects endogenous levels of total TNNI2
Immunogen Description	A synthesized peptide derived from human TNNI2
Other Names	DA2B; FSSV; fsTnl; AMCD2B;
Accession No.	Uniprot:P48788
Uniprot	P48788
Calculated MW	21kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:1000~1:5000IP:1:30

Images



Western blot analysis of Troponin I fast skeletal muscle expression in Human skeletal muscle lysate.

Product Description

Troponin I is the inhibitory subunit of troponin, the thin filament regulatory complex which confers calcium-sensitivity to striated muscle actomyosin ATPase activity.

Note: This product is for in vitro research use only