

ACTN3 Rabbit mAb

Catalog No: #56519

Package Size: #56519-1 50ul #56519-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ACTN3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	ACTN3 Antibody detects endogenous levels of total ACTN3
Immunogen Description	A synthesized peptide derived from human ACTN3
Other Names	Actinin alpha 3; ACTN3;
Accession No.	Uniprot:Q08043
Uniprot	Q08043
Calculated MW	103 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of ACTN3 expression in mouse skeletal muscle lysate.

Product Description

F-actin cross-linking protein which is thought to anchor actin to a variety of intracellular structures. This is a bundling protein.

Note: This product is for in vitro research use only