

RAIDD Rabbit mAb

Catalog No: #56520

Package Size: #56520-1 50ul #56520-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	RAIDD Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC IP FC
Species Reactivity	Human Rat
Specificity	RAIDD Antibody detects endogenous levels of total RAIDD
Immunogen Description	A synthesized peptide derived from human RAIDD
Other Names	Cradd;
Accession No.	Uniprot:P78560
Uniprot	P78560
Calculated MW	23 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details


WB: 1:500~1:2000

ICC/IF: 1:50~1:200

IP: 1:50

FC: 1:50

Images



Western blot analysis of RAIDD expression in HeLa cell lysate.

Product Description

Apoptotic adaptor molecule specific for caspase-2 and FASL/TNF receptor-interacting protein RIP. In the presence of RIP and TRADD, CRADD

recruits caspase-2 to the TNFR-1 signalling complex.

Note: This product is for in vitro research use only