

TPPP Rabbit mAb

Catalog No: #56524

Package Size: #56524-1 50ul #56524-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TPPP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	TPPP Antibody detects endogenous levels of total TPPP
Immunogen Description	A synthesized peptide derived from human TPPP
Other Names	p24; p25; p25alpha; TPPP; TPPP/p25; TPPP1; Tubulin polymerization promoting protein;
Accession No.	Uniprot:O94811
Uniprot	O94811
Calculated MW	25 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200FC:1:50

Images



Western blot analysis of TPPP expression in SHSY5Y cell lysate.

Product Description

May play a role in the polymerization of tubulin into microtubules, microtubule bundling and the stabilization of existing microtubules, thus maintaining the integrity of the microtubule network. May play a role in mitotic spindle assembly and nuclear envelope breakdown.

Note: This product is for in vitro research use only