

PON1 Rabbit mAb

Catalog No: #56538

Package Size: #56538-1 50ul #56538-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PON1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human
Specificity	PON1 Antibody detects endogenous levels of total PON1
Immunogen Description	A synthesized peptide derived from human PON1
Other Names	A esterase 1; Arylesterase 1; ESA; MVCD5; Paraoxonase 1; PON1;
Accession No.	Uniprot:P27169
Uniprot	P27169
Calculated MW	40 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200IP:1:50

Images



Western blot analysis of PON1 expression in human plasma cell lysate.

Product Description

Hydrolyzes the toxic metabolites of a variety of organophosphorus insecticides. Capable of hydrolyzing a broad spectrum of organophosphate substrates and lactones, and a number of aromatic carboxylic acid esters.

Note: This product is for in vitro research use only