

PCYT1A Rabbit mAb

Catalog No: #56541

Package Size: #56541-1 50ul #56541-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PCYT1A Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human
Specificity	PCYT1A Antibody detects endogenous levels of total PCYT1A
Immunogen Description	A synthesized peptide derived from human PCYT1A
Other Names	CCT A; CCT alph; CTPCT; PCYT1; Pcyt1a;
Accession No.	Uniprot:P49585
Uniprot	P49585
Calculated MW	42 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000 ICC/IF: 1:50~1:200 IP: 1:50 FC: 1:50

Images



Western blot analysis of PCYT1A expression in A549 cell lysate.

Product Description

Controls phosphatidylcholine synthesis.

Note: This product is for in vitro research use only