

Cripto1 Rabbit mAb

Catalog No: #56545

Package Size: #56545-1 50ul #56545-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Cripto1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human
Specificity	Cripto1 Antibody detects endogenous levels of total Cripto1
Immunogen Description	A synthesized peptide derived from human Cripto1
Other Names	CRGF; Cripto 1; cripto; Cripto1 growth factor; TDGF1;
Accession No.	Uniprot:P13385
Uniprot	P13385
Calculated MW	21 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IC/IF:1:50~1:200IP:1:50FC:1:50

Images



Western blot analysis of Cripto1 expression in NCCIT cell lysate.

Product Description

Could play a role in the determination of the epiblastic cells that subsequently give rise to the mesoderm.

Note: This product is for in vitro research use only