

## APRIL Rabbit mAb

Catalog No: #56549

Package Size: #56549-1 50ul #56549-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	APRIL Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC
Species Reactivity	Human Mouse Rat
Specificity	APRIL Antibody detects endogenous levels of total APRIL
Immunogen Description	A synthesized peptide derived from human APRIL
Other Names	APRIL; CD256; TALL2; TNFSF13; TRDL1; ZTNF2;
Accession No.	Uniprot:O75888
Uniprot	O75888
Calculated MW	35kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB: 1:500~1:2000 ICC/IF: 1:50~1:200

## Images



Western blot analysis of APRIL expression in Jurkat cell lysate.

## Product Description

Cytokine that binds to TNFRSF13B/TACI and to TNFRSF17/BCMA. May be implicated in the regulation of tumor cell growth. May be involved in monocyte/macrophage-mediated immunological processes.

---

Note: This product is for in vitro research use only