

CD97 Rabbit mAb

Catalog No: #56568

Package Size: #56568-1 50ul #56568-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CD97 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC FC
Species Reactivity	Human Rat
Specificity	CD97 Antibody detects endogenous levels of total CD97
Immunogen Description	A synthesized peptide derived from human CD97
Other Names	ADGRE5; CD97; TM7LN1;
Accession No.	Uniprot:P48960
Uniprot	P48960
Calculated MW	55/28kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.


Application Details

WB: 1:500~1:2000

IHC: 1:50~1:200


FC: 1:50

Images



Western blot analysis of CD97 expression in human tonsil cell lysate.

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Product Description

Receptor potentially involved in both adhesion and signaling processes early after leukocyte activation. Plays an essential role in leukocyte migration.

Note: This product is for in vitro research use only