

## NeuroD2 Rabbit mAb

Catalog No: #56569

Package Size: #56569-1 50ul #56569-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	NeuroD2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	NeuroD2 Antibody detects endogenous levels of total NeuroD2
Immunogen Description	A synthesized peptide derived from human NeuroD2
Other Names	bHLHa1; NDR2; NDRF; neuroD related factor; NeuroD2; neurogenic differentiation 2;
Accession No.	Uniprot:Q15784
Uniprot	Q15784
Calculated MW	48kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB: 1:500~1:2000

## Images



Western blot analysis of NeuroD2 expression in mouse cerebellum cell lysate.

## Product Description

Appears to mediate neuronal differentiation.

---

Note: This product is for in vitro research use only