

BMP15 Rabbit mAb

Catalog No: #56576

Package Size: #56576-1 50ul #56576-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	BMP15 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human Mouse Rat
Specificity	BMP15 Antibody detects endogenous levels of total BMP15
Immunogen Description	A synthesized peptide derived from human BMP15
Other Names	BMP15; GDF9B; Growth/differentiation factor 9B; ODG2; POF4;
Accession No.	Uniprot:O95972
Uniprot	O95972
Calculated MW	50kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000FC:1:50

Images



Western blot analysis of BMP15 expression in HeLa cell lysate.

Product Description

May be involved in follicular development. Oocyte-specific growth/differentiation factor that stimulates folliculogenesis and granulosa cell (GC) growth.

Note: This product is for in vitro research use only