

TIMM50 Rabbit mAb

Catalog No: #56579

Package Size: #56579-1 50ul #56579-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TIMM50 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse
Specificity	TIMM50 Antibody detects endogenous levels of total TIMM50
Immunogen Description	A synthesized peptide derived from human TIMM50
Other Names	TIM50; TIM50L; TIMM50;
Accession No.	Uniprot:Q3ZCQ8
Uniprot	Q3ZCQ8
Calculated MW	36kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000 IHC: 1:50~1:200

Images



Western blot analysis of TIMM50 expression in 293T cell lysate.

Product Description

Essential component of the TIM23 complex, a complex that mediates the translocation of transit peptide-containing proteins across the mitochondrial inner membrane. Has some phosphatase activity in vitro; however such activity may not be relevant in vivo.

Note: This product is for in vitro research use only