

BCL2L12 Rabbit mAb

Catalog No: #56580

Package Size: #56580-1 50ul #56580-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	BCL2L12 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human
Specificity	BCL2L12 Antibody detects endogenous levels of total BCL2L12
Immunogen Description	A synthesized peptide derived from human BCL2L12
Other Names	BCL2 like 12; BCL2L12; BPR;
Accession No.	Uniprot:Q9HB09
Uniprot	Q9HB09
Calculated MW	32kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000


IHC: 1:50~1:200

ICC/IF: 1:50~1:200

IP: 1:50


FC: 1:50

Images



Western blot analysis of BCL2L12 expression in MCF7 cell lysate.

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Product Description

Belongs to the Bcl-2 family.

Note: This product is for in vitro research use only