

NEK6 Rabbit mAb

Catalog No: #56581

Package Size: #56581-1 50ul #56581-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NEK6 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human Mouse Rat
Specificity	NEK6 Antibody detects endogenous levels of total NEK6
Immunogen Description	A synthesized peptide derived from human NEK6
Other Names	NEK6; NimA related protein kinase 6; SID61512;
Accession No.	Uniprot:Q9HC98
Uniprot	Q9HC98
Calculated MW	28kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000 IP: 1:50

Images



Western blot analysis of NEK6 expression in HeLa cell lysate.

Product Description

Activated during M phase. Required for chromosome segregation at metaphase-anaphase transition and therefore for mitotic progression. Inhibition of activity results in apoptosis. Phosphorylates KIF11 to promote mitotic spindle formation.

Note: This product is for in vitro research use only