

LTA4H Rabbit mAb

Catalog No: #56583

Package Size: #56583-1 50ul #56583-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	LTA4H Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	LTA4H Antibody detects endogenous levels of total LTA4H
Immunogen Description	A synthesized peptide derived from human LTA4H
Other Names	LTA4; LTA4 hydrolase; LTA4H;
Accession No.	Uniprot:P09960
Uniprot	P09960
Calculated MW	69kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200FC:1:50

Images



Western blot analysis of LTA4H expression in A549 cell lysate.

Product Description

Hydrolyzes an epoxide moiety of leukotriene A4 (LTA-4) to form leukotriene B4 (LTB-4). The enzyme also has some peptidase activity.

Note: This product is for in vitro research use only