

## FANCA Rabbit mAb

Catalog No: #56588

Package Size: #56588-1 50ul #56588-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	FANCA Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human
Specificity	FANCA Antibody detects endogenous levels of total FANCA
Immunogen Description	A synthesized peptide derived from human FANCA
Other Names	FA 1; FA; FA H; FA1; FAA; FACA; FAH; Fanca; FANCH;
Accession No.	Uniprot:O15360
Uniprot	O15360
Calculated MW	163Da
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:1000IP:1:40

## Images



Western blot analysis of FANCA expression in A431 cell lysate.

## Product Description

DNA repair protein that may operate in a postreplication repair or a cell cycle checkpoint function. May be involved in interstrand DNA cross-link repair and in the maintenance of normal chromosome stability.

---

Note: This product is for in vitro research use only