

RGS18 Rabbit mAb

Catalog No: #56594

Package Size: #56594-1 50ul #56594-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	RGS18 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	RGS18 Antibody detects endogenous levels of total RGS18
Immunogen Description	A synthesized peptide derived from human RGS18
Other Names	RGS13; RGS18;
Accession No.	Uniprot:Q9NS28
Uniprot	Q9NS28
Calculated MW	28kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of RGS18 expression in Human PBMC lysate.

Product Description

Inhibits signal transduction by increasing the GTPase activity of G protein alpha subunits thereby driving them into their inactive GDP-bound form.

Binds to G(i) alpha-1, G(i) alpha-2, G(i) alpha-3 and G(q) alpha.

Note: This product is for in vitro research use only