

Oxytocin Receptor Rabbit mAb

Catalog No: #56597

Package Size: #56597-1 50ul #56597-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Oxytocin Receptor Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	Oxytocin Receptor Antibody detects endogenous levels of total Oxytocin Receptor
Immunogen Description	A synthesized peptide derived from human Oxytocin Receptor
Other Names	OT R; OTR; OTR1; OXT R; OxtR;
Accession No.	Uniprot:P30559
Uniprot	P30559
Calculated MW	43kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of Oxytocin Receptor expression in (1) Jurkat cell lysate; (2) Mouse spleen lysate.

Product Description

Receptor for oxytocin. The activity of this receptor is mediated by G proteins which activate a phosphatidylinositol-calcium second messenger system.

Note: This product is for in vitro research use only