

## ZA2G / AZGP1 Rabbit mAb

Catalog No: #56604

Package Size: #56604-1 50ul #56604-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	ZA2G / AZGP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human
Specificity	ZA2G / AZGP1 Antibody detects endogenous levels of total ZA2G / AZGP1
Immunogen Description	A synthesized peptide derived from human ZA2G / AZGP1
Other Names	Azgp1; ZA2G; ZAG; ZNGP1;
Accession No.	Uniprot:P25311
Uniprot	P25311
Calculated MW	43kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IHC:1:50~1:200IP:1:40

## Images



Western blot analysis of ZA2G / AZGP1 expression in HepG2 cell lysate.

## Product Description

Stimulates lipid degradation in adipocytes and causes the extensive fat losses associated with some advanced cancers. May bind polyunsaturated fatty acids.

---

Note: This product is for in vitro research use only