

## Neuronatin Rabbit mAb

Catalog No: #56613

Package Size: #56613-1 50ul #56613-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Neuronatin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse
Specificity	Neuronatin Antibody detects endogenous levels of total Neuronatin
Immunogen Description	A synthesized peptide derived from human Neuronatin
Other Names	Neuronatin; NNAT; Peg 5; Peg5;
Accession No.	Uniprot:Q16517
Uniprot	Q16517
Calculated MW	9kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB: 1:500~1:2000

## Images



Western blot analysis of Neuronatin expression in SH-SY5Y cell lysate.

## Product Description

May participate in the maintenance of segment identity in the hindbrain and pituitary development, and maturation or maintenance of the overall structure of the nervous system. May function as a regulatory subunit of ion channels.

---

Note: This product is for in vitro research use only