

VPS18 Rabbit mAb

Catalog No: #56621

Package Size: #56621-1 50ul #56621-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	VPS18 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	VPS18 Antibody detects endogenous levels of total VPS18
Immunogen Description	A synthesized peptide derived from human VPS18
Other Names	hVPS18; PEP3; vps18;
Accession No.	Uniprot:Q9P253
Uniprot	Q9P253
Calculated MW	100kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000ICC/IF:1:50~1:200FC:1:50

Images



Western blot analysis of VPS18 expression in K562 cell lysate.

Product Description

May play a role in vesicle-mediated protein trafficking to lysosomal compartments and in membrane docking/fusion reactions of late endosomes/lysosomes.

Note: This product is for in vitro research use only