

CCL19 Rabbit mAb

Catalog No: #56633

Package Size: #56633-1 50ul #56633-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CCL19 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human
Specificity	CCL19 Antibody detects endogenous levels of total CCL19
Immunogen Description	A synthesized peptide derived from human CCL19
Other Names	CCL19; CK beta 11; CKb11; ELC; Exodus3; MIP3 beta; MIP3B; SCYA19;
Accession No.	Uniprot:Q99731
Uniprot	Q99731
Calculated MW	11kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:1000~1:5000FC:1:100

Images



Western blot analysis of CCL19 expression in Human Macrophage Inflammatory Protein 3 beta recombinant protein lysate.

Product Description

May play a role not only in inflammatory and immunological responses but also in normal lymphocyte recirculation and homing. May play an important role in trafficking of T-cells in thymus, and T-cell and B-cell migration to secondary lymphoid organs.

Note: This product is for in vitro research use only