

PER2 Rabbit mAb

Catalog No: #56641

Package Size: #56641-1 50ul #56641-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PER2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human
Specificity	PER2 Antibody detects endogenous levels of total PER2
Immunogen Description	A synthesized peptide derived from human PER2
Other Names	FASPS; FASPS1; hPER2; hPER2; KIAA0347; PER2; Period 2;
Accession No.	Uniprot:O15055
Uniprot	O15055
Calculated MW	137kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000ICC/IF:1:50~1:200FC:1:50

Images



Western blot analysis of PER2 expression in HeLa cell lysate.

Product Description

Component of the circadian clock mechanism which is essential for generating circadian rhythms. Negative element in the circadian transcriptional loop. Influences clock function by interacting with other circadian regulatory proteins and transporting them to the nucleus. Negatively regulates CLOCK

Note: This product is for in vitro research use only