

DDX17 Rabbit mAb

Catalog No: #56642

Package Size: #56642-1 50ul #56642-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	DDX17 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse
Specificity	DDX17 Antibody detects endogenous levels of total DDX17
Immunogen Description	A synthesized peptide derived from human DDX17
Other Names	Ddx17; P72; RH70;
Accession No.	Uniprot:Q92841
Uniprot	Q92841
Calculated MW	70,80kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200FC:1:100

Images



Western blot analysis of DDX17 expression in 293T cell lysate.

Product Description

RNA-dependent ATPase activity.

Note: This product is for in vitro research use only