

CRYAA Rabbit mAb

Catalog No: #56645

Package Size: #56645-1 50ul #56645-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CRYAA Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	CRYAA Antibody detects endogenous levels of total CRYAA
Immunogen Description	A synthesized peptide derived from human CRYAA
Other Names	Acry 1; CRYA1; CRYAA; HspB4;
Accession No.	Uniprot:P02489
Uniprot	P02489
Calculated MW	17kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000ICC/IF:1:50~1:200

Images



Western blot analysis of CRYAA expression in Mouse eyeball lysate.

Product Description

May contribute to the transparency and refractive index of the lens.

Note: This product is for in vitro research use only