

## MDH1 Rabbit mAb

Catalog No: #56648

Package Size: #56648-1 50ul #56648-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	MDH1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	MDH1 Antibody detects endogenous levels of total MDH1
Immunogen Description	A synthesized peptide derived from human MDH1
Other Names	MDH s; mdh1; MDHA; MDHs; MOR2;
Accession No.	Uniprot:P40925
Uniprot	P40925
Calculated MW	36kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:1000~1:5000|ICC/IF:1:50~1:200|IP:1:50|FC:1:30

## Images



Western blot analysis of MDH1 expression in HeLa cell lysate.

## Product Description

Belongs to the LDH/MDH superfamily. MDH type 2 family.

---

Note: This product is for in vitro research use only