

Granulin Rabbit mAb

Catalog No: #56657

Package Size: #56657-1 50ul #56657-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Granulin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human
Specificity	Granulin Antibody detects endogenous levels of total Granulin
Immunogen Description	A synthesized peptide derived from human Granulin
Other Names	Acrogranin; CLN11; GEP; GP88; Granulins; GRN; PCDGF; PEPI; PGRN; Proepithelin; Progranulin;
Accession No.	Uniprot:P28799
Uniprot	P28799
Calculated MW	55,75kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:50FC:1:100

Images



Western blot analysis of Granulin expression in 293T cell lysate.

Product Description

Granulins have possible cytokine-like activity. They may play a role in inflammation, wound repair, and tissue remodeling.

Note: This product is for in vitro research use only