

GPCR / LGR6 Rabbit mAb

Catalog No: #56659

Package Size: #56659-1 50ul #56659-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GPCR / LGR6 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	GPCR / LGR6 Antibody detects endogenous levels of total GPCR / LGR6
Immunogen Description	A synthesized peptide derived from human GPCR / LGR6
Other Names	GPCR; Lgr6;
Accession No.	Uniprot:Q9HBX8
Uniprot	Q9HBX8
Calculated MW	104kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200FC:1:50

Images



Western blot analysis of GPCR / LGR6 expression in HUVEC cell lysate.

Product Description

Orphan receptor.

Note: This product is for in vitro research use only