

# ICOS Rabbit mAb

Catalog No: #56660

Package Size: #56660-1 50ul #56660-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ICOS Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	ICOS Antibody detects endogenous levels of total ICOS
Immunogen Description	A synthesized peptide derived from human ICOS
Other Names	AILIM; CD278; CRP1; CVID1; ICOS;
Accession No.	Uniprot:Q9Y6W8
Uniprot	Q9Y6W8
Calculated MW	22kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000ICC/IF:1:50~1:200FC:1:50

## Images



Western blot analysis of ICOS expression in Jurkat cell lysate.

## Product Description

Enhances all basic T-cell responses to a foreign antigen, namely proliferation, secretion of lymphokines, up-regulation of molecules that mediate cell-cell interaction, and effective help for antibody secretion by B-cells.

---

Note: This product is for in vitro research use only