

LDHB Rabbit mAb

Catalog No: #56661

Package Size: #56661-1 50ul #56661-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|--|
| Product Name | LDHB Rabbit mAb |
| Host Species | Rabbit |
| Clonality | Monoclonal |
| Isotype | Rabbit IgG |
| Purification | Affinity-chromatography |
| Applications | WB IHC ICC/IF IP FC |
| Species Reactivity | Human |
| Specificity | LDHB Antibody detects endogenous levels of total LDHB |
| Immunogen Description | A synthesized peptide derived from human LDHB |
| Other Names | HEL S 281; LDH B; LDH H; LDH-B; LDH-H; LDHB; LDHBD; LDHH; TRG5; |
| Accession No. | Uniprot:P07195 |
| Uniprot | P07195 |
| Calculated MW | 37kDa |
| Formulation | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage | Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle. |

Application Details

WB:1:1000~1:5000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:100FC:1:60

Images



Western blot analysis of LDHB expression in HeLa cell lysate.

Product Description

Fermentation; pyruvate fermentation to lactate; (S)-lactate from pyruvate: step 1/1.

Note: This product is for in vitro research use only