

HOXA13 Rabbit mAb

Catalog No: #56671

Package Size: #56671-1 50ul #56671-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	HOXA13 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	HOXA13 Antibody detects endogenous levels of total HOXA13
Immunogen Description	A synthesized peptide derived from human HOXA13
Other Names	Homeobox1J; HomeoboxA13; HOX 1; HOX 1J; HOX A13; HOX1; HOX1J; HOXA 13; HOXA13; Transcription factor HOXA13;
Accession No.	Uniprot:P31271
Uniprot	P31271
Calculated MW	42kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

Images



Western blot analysis of HOXA13 expression in 293T cell lysate.

Product Description

Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.

Note: This product is for in vitro research use only