

VPS24 Rabbit mAb

Catalog No: #56672

Package Size: #56672-1 50ul #56672-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	VPS24 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP FC
Species Reactivity	Human Mouse Rat
Specificity	VPS24 Antibody detects endogenous levels of total VPS24
Immunogen Description	A synthesized peptide derived from human VPS24
Other Names	VPS24; CGI149; CHMP3; hVps24; NEDF;
Accession No.	Uniprot:Q9Y3E7
Uniprot	Q9Y3E7
Calculated MW	29kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200IP:1:50FC:1:50

Images



Western blot analysis of VPS24 expression in HepG2 cell lysate.

Product Description

Probable core component of the endosomal sorting required for transport complex III (ESCRT-III) which is involved in multivesicular bodies (MVBs) formation and sorting of endosomal cargo proteins into MVBs.

Note: This product is for in vitro research use only