

NCF1 Rabbit mAb

Catalog No: #56681

Package Size: #56681-1 50ul #56681-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NCF1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human
Specificity	NCF1 Antibody detects endogenous levels of total NCF1
Immunogen Description	A synthesized peptide derived from human NCF1
Other Names	NCF1; NCF1A; NOXO2; SH3PXD1A; p47phox;
Accession No.	Uniprot:P14598
Uniprot	P14598
Calculated MW	45kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000 IP: 1:50

Images



Western blot analysis of NCF1 expression in Raji cell lysate.

Product Description

NCF2, NCF1, and a membrane bound cytochrome b558 are required for activation of the latent NADPH oxidase (necessary for superoxide production).

Note: This product is for in vitro research use only