

## PDZK1 Rabbit mAb

Catalog No: #56686

Package Size: #56686-1 50ul #56686-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	PDZK1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human
Specificity	PDZK1 Antibody detects endogenous levels of total PDZK1
Immunogen Description	A synthesized peptide derived from human PDZK1
Other Names	CAP70; CLAMP; PDZD1; NHERF3; NHERF-3;
Accession No.	Uniprot:Q5T2W1
Uniprot	Q5T2W1
Calculated MW	70kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IP:1:50

## Images



Western blot analysis of PDZK1 expression in T47-D cell lysate.

## Product Description

A scaffold protein that connects plasma membrane proteins and regulatory components, regulating their surface expression in epithelial cells apical domains.

---

Note: This product is for in vitro research use only