

## Pumilio 2 Rabbit mAb

Catalog No: #56688

Package Size: #56688-1 50ul #56688-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Pumilio 2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	Pumilio 2 Antibody detects endogenous levels of total Pumilio 2
Immunogen Description	A synthesized peptide derived from human Pumilio 2
Other Names	PUM2; PUMH2; Pumilio2; PUML2;
Accession No.	Uniprot:Q8TB72
Uniprot	Q8TB72
Calculated MW	114kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000ICC/IF:1:50~1:200IP:1:20FC:1:50

## Images



Western blot analysis of Pumilio 2 expression in 293T cell lysate.

## Product Description

Sequence-specific RNA-binding protein that regulates translation and mRNA stability by binding the 3'-UTR of mRNA targets. Its interactions and tissue specificity suggest that it may be required to support proliferation and self-renewal of stem cells by regulating the translation of key transcripts.

---

Note: This product is for in vitro research use only