

Engrailed 1 Rabbit mAb

Catalog No: #56692

Package Size: #56692-1 50ul #56692-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Engrailed 1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP FC
Species Reactivity	Human Mouse Rat
Specificity	Engrailed 1 Antibody detects endogenous levels of total Engrailed 1
Immunogen Description	A synthesized peptide derived from human Engrailed 1
Other Names	EN1; Engrailed 1; Engrailed; Hu-En-1; Mo-en.1;
Accession No.	Uniprot:Q05925
Uniprot	Q05925
Calculated MW	55kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IP:1:50FC:1:100

Images



Western blot analysis of Engrailed 1 expression in 293T cell lysate.

Product Description

Belongs to the engrailed homeobox family.

Note: This product is for in vitro research use only