

TF2B Rabbit mAb

Catalog No: #56697

Package Size: #56697-1 50ul #56697-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TF2B Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human Mouse Rat
Specificity	TF2B Antibody detects endogenous levels of total TF2B
Immunogen Description	A synthesized peptide derived from human TF2B
Other Names	gtf2b; S300 II; TF IIB; TF2B; TFIIB;
Accession No.	Uniprot:Q00403
Uniprot	Q00403
Calculated MW	31kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IP:1:50

Images



Western blot analysis of TFIIB expression in (1) HepG2 cell lysate; (2) RAW264.7 cell lysate.

Product Description

General factor that plays a major role in the activation of eukaryotic genes transcribed by RNA polymerase II.

Note: This product is for in vitro research use only