

PAR4 Rabbit mAb

Catalog No: #56698

Package Size: #56698-1 50ul #56698-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PAR4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	PAR4 Antibody detects endogenous levels of total PAR4
Immunogen Description	A synthesized peptide derived from human PAR4
Other Names	PAR 4; Pawr;
Accession No.	Uniprot:Q96IZ0
Uniprot	Q96IZ0
Calculated MW	40kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of PAR4 expression in LnCaP cell lysate.

Product Description

Pro-apoptotic protein capable of selectively inducing apoptosis in cancer cells, sensitizing the cells to diverse apoptotic stimuli and causing regression of tumors in animal models. Induces apoptosis in certain cancer cells by activation of the Fas prodeath pathway and coparallel inhibition of NF-kappa-B transcriptional activity.

Note: This product is for in vitro research use only