

ADA Rabbit mAb

Catalog No: #56699

Package Size: #56699-1 50ul #56699-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ADA Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP
Species Reactivity	Human Rat
Specificity	ADA Antibody detects endogenous levels of total ADA
Immunogen Description	A synthesized peptide derived from human ADA
Other Names	ada; ADA1;
Accession No.	Uniprot:P00813
Uniprot	P00813
Calculated MW	41kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000ICC/IF:1:50~1:200IP:1:50

Images



Western blot analysis of ADA expression in Jurkat cell lysate.

Product Description

Catalyzes the hydrolytic deamination of adenosine and 2-deoxyadenosine. Plays an important role in purine metabolism and in adenosine homeostasis. Modulates signaling by extracellular adenosine, and so contributes indirectly to cellular signaling events.

Note: This product is for in vitro research use only