

MALT1 Rabbit mAb

Catalog No: #56700

Package Size: #56700-1 50ul #56700-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	MALT1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human
Specificity	MALT1 Antibody detects endogenous levels of total MALT1
Immunogen Description	A synthesized peptide derived from human MALT1
Other Names	Malt1; MLT; MLT1; Paracaspase;
Accession No.	Uniprot:Q9UDY8
Uniprot	Q9UDY8
Calculated MW	92kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

ICC/IF: 1:50~1:200

Images



Western blot analysis of MALT1 expression in K562 cell lysate.

Product Description

Enhances BCL10-induced activation of NF-kappa-B. Involved in nuclear export of BCL10. Binds to TRAF6, inducing TRAF6 oligomerization and activation of its ligase activity.

Note: This product is for in vitro research use only