

HLA B7 Rabbit mAb

Catalog No: #56703

Package Size: #56703-1 50ul #56703-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	HLA B7 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	HLA B7 Antibody detects endogenous levels of total HLA B7
Immunogen Description	A synthesized peptide derived from human HLA B7
Other Names	HLA B; HLAB;
Accession No.	Uniprot:P01889
Uniprot	P01889
Calculated MW	40kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200

Images



Western blot analysis of HLA B7 expression in Ramos cell lysate.

Product Description

HLA B belongs to the HLA class I heavy chain paralogues. This class I molecule is a heterodimer consisting of a heavy chain and a light chain (beta 2 microglobulin). The heavy chain is anchored in the membrane.

Note: This product is for in vitro research use only