

## MYF5 Rabbit mAb

Catalog No: #56708

Package Size: #56708-1 50ul #56708-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	MYF5 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human Mouse Rat
Specificity	MYF5 Antibody detects endogenous levels of total MYF5
Immunogen Description	A synthesized peptide derived from human MYF5
Other Names	Myf5; Class C basic helix loop helix protein 2
Accession No.	Uniprot:P13349
Uniprot	P13349
Calculated MW	37kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:1000~1:5000IP:1:50

## Images



Western blot analysis of MYF5 expression in Human skeletal muscle lysate.

## Product Description

Involved in muscle differentiation (myogenic factor). Induces fibroblasts to differentiate into myoblasts. Probable sequence specific DNA-binding protein.

---

Note: This product is for in vitro research use only