

## S1PR3 Rabbit mAb

Catalog No: #56710

Package Size: #56710-1 50ul #56710-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	S1PR3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	S1PR3 Antibody detects endogenous levels of total S1PR3
Immunogen Description	A synthesized peptide derived from human S1PR3
Other Names	EDG3; LPB3; S1P3; S1pr3;
Accession No.	Uniprot:Q99500
Uniprot	Q99500
Calculated MW	42kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle.

## Application Details

WB: 1:500~1:2000

## Images



Western blot analysis of S1PR3 expression in HepG2 cell lysate.

## Product Description

Receptor for the lysosphingolipid sphingosine 1-phosphate (S1P). S1P is a bioactive lysophospholipid that elicits diverse physiological effect on most types of cells and tissues.

---

Note: This product is for in vitro research use only