

Thyroid Peroxidase Rabbit mAb

Catalog No: #56712

Package Size: #56712-1 50ul #56712-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Thyroid Peroxidase Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human
Specificity	Thyroid Peroxidase Antibody detects endogenous levels of total Thyroid Peroxidase
Immunogen Description	A synthesized peptide derived from human Thyroid Peroxidase
Other Names	MSA; TDH2A; TPO; TPX;
Accession No.	Uniprot:P07202
Uniprot	P07202
Calculated MW	103kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:1000~1:5000IHC:1:100~1:500IP:1:50

Images



Western blot analysis of Thyroid Peroxidase expression in Human thyroid lysate.

Product Description

Iodination and coupling of the hormonogenic tyrosines in thyroglobulin to yield the thyroid hormones T(3) and T(4).

Note: This product is for in vitro research use only