

## IL17A Receptor Rabbit mAb

Catalog No: #56715

Package Size: #56715-1 50ul #56715-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	IL17A Receptor Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human
Specificity	IL17A Receptor Antibody detects endogenous levels of total IL17A Receptor
Immunogen Description	A synthesized peptide derived from human IL17A Receptor
Other Names	IL-17 receptor; IL-17RA; IL17RA; interleukin-17 receptor;
Accession No.	Uniprot:Q96F46
Uniprot	Q96F46
Calculated MW	160kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:1000~1:5000IP:1:50

## Images



Western blot analysis of IL17A Receptor expression in Raji cell lysate.

## Product Description

Receptor for IL17A and IL17F. Binds its IL17A ligand with low affinity, suggesting that additional components are involved in IL17A-induced signaling.

---

Note: This product is for in vitro research use only