

CHRNA5 Rabbit mAb

Catalog No: #56723

Package Size: #56723-1 50ul #56723-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CHRNA5 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP FC
Species Reactivity	Human Mouse Rat
Specificity	CHRNA5 Antibody detects endogenous levels of total CHRNA5
Immunogen Description	A synthesized peptide derived from human CHRNA5
Other Names	AChR; CHRNA5; LNCR2; NACHRA5;
Accession No.	Uniprot:P30532
Uniprot	P30532
Calculated MW	55kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IP:1:50FC:1:50

Images



Western blot analysis of CHRNA5 expression in A431 cell lysate.

Product Description

After binding acetylcholine, the AChR responds by an extensive change in conformation that affects all subunits and leads to opening of an ion-conducting channel across the plasma membrane.

Note: This product is for in vitro research use only